
ECO II+

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1. Notes

1.1 Translations

1	Translation of the original instructions for use
---	--

1.2 Images

The images contained in these instructions for use are for reference and may differ from the actual appearance of the product.

1.3 General instructions

METASYS can only guarantee the safety, reliability and performance of the dental device if the following instructions are adhered to:

- The product is only to be used in accordance with the instructions for use.
- During maintenance and service work (inspection, service, repair, replacement), only use original replacement parts.
- All manufacturer specifications for the treatment units to which the device is connected must be observed.
- After commissioning, complete the proof of installation and send this to METASYS in order to define the warranty period.
- All maintenance and service work must be entered into the device logbook.
- On request by an authorised technician, METASYS will provide all documentation that may be of use to technically qualified personnel during maintenance and service works
- METASYS accepts no responsibility for damage that may arise due to external influences (defective installation), using incorrect information, improper use of the dental device, or maintenance and service works being carried out improperly.
- The user must familiarise himself/herself with how to operate the dental device and ensure that the dental device is in good condition each time before operating it.
- Important: Read the accompanying documents of the device carefully before installation, commissioning and use and keep them for the entire service life of the product!

1.4 Explanation of the symbols

2	Item number	13	Use eye protection
3	Serial number	14	Use mouth and nose protection
4	Manufacturer	15	Information
5	Date of manufacture	16	Protect from heat / Protect from sunlight
6	Follow instructions for use	17	Protect from moisture / Store in a dry place
7	Temperature limit	18	Fragile
8	Humidity limit	19	Package orientation at the top
9	Maximum installation altitude	20	Mark of conformity per para. 22 section 4 Lower Saxony Building Regs
10	Caution / attention		
11	Responsible for the translation		
12	Use hand protection		

1.5 Safety-related information

All serious incidents related to the device must be reported to the manufacturer and the competent authority of the Member State where the user and/or the patient is resident.

2. Product description

The METASYS ECO II+ is a sedimentation amalgam separator for the installation after a central suction system or in the waste water pipe of the dental unit. When used in a liquid ring system, the ECO II+ can be installed in the suction line before the liquid ring pump.

2.1 Technical data

Max. ambient temperature	40 °C
Separation rate	≥ 95%
Max. water flow rate (via air/water separation)	1 l / min
Collection container volume	700 cm ³
Empty weight ECO II+ with holder	1,7 kg 3,2 kg*
Full weight ECO II+ container	2,7 kg 4,2 kg*
Weight of cover (without ECO II+)	8 kg
21 - Dimensions (H x W x D)	309 x 163 x 190 mm
Dimensions with cover (H x W x D)	342 x 218 x 260 mm
DIBT approval number	Z-64.3-26

Possible suction systems:

- dry-line vacuum generators with air/water separation, e.g. central suction system METASYS EX-COM hybrid.
- liquid ring systems

With the ECO II+ a separation rate of 95% is achieved. The use of air abrasion granulate does not affect the separation rate of the amalgam separator!

2.2 Type plate

22 - The type plate is located on module 1.

A	Device description	F	Date of manufacture
B	Item number	G	Manufacturer
C	Serial number	H	Mark of conformity per para. 22 section 4 Lower Saxony Building Regs
D	Connection data	I	Follow instructions for use
E	DIBT approval number		

2.3 Construction

22 The amalgam separator ECO II+ consists of 2 modules:

22-1 Module 1 is the central mounting element and at the same time wall bracket.

22-2 Module 2 is the sedimentation amalgam separator with water connections.

2.4 Functional description

The amalgam separator ECO II+ works on the principle of sedimentation, divided into several sections. The waste water from the dental units runs from the central suction system with air/water separation into the inlet nozzle of the ECO II+. In liquid ring systems, the water flow from the suction line into the inlet nozzle. From there the waste water is led through the first large sedimentation chamber. Afterwards it is led through the first sedimentation cartridge, through cascades and the second sedimentation cartridge into the second sedimentation chamber. Through the third sedimentation cartridge the cleaned water is led directly into the sewage pipe. In liquid ring systems, the clean water is led into the liquid ring pump.

23

A	inlet nozzle	E	second sedimentation cartridge
B	first sedimentation chamber	F	second sedimentation chamber
C	first sedimentation cartridge	G	third sedimentation cartridge
D	cascade	H	waste water connection

Product appearance:

Please note that you may see small loose silver-colored particles (zinc) inside the ECO II+ amalgam separator. These particles do not affect the cartridges or performance of the ECO II+. If you have any questions, please feel free to contact METASYS Medizintechnik GmbH.

3. Preparation for use

3.1 Transport and storage

ECO II+ is shipped in a cardboard box with two polystyrene trays. This original packaging must be used for any onward and return transport. The device must always be transported and stored in an upright position. The device must be transported to the installation site in a completely packed state. After unpacking the device, check for completeness and possible transport damage.

24	Transport and storage	25	Installation requirements
-----------	-----------------------	-----------	---------------------------

3.2 Installation options

The ECO II+ allows a variety of different installation options. It is necessary to check if the possible water input from simultaneously operated units can exceed 1 l/min.

If the maximal amount of water can exceed 1 l/min or if the waste water from several units is brought together, then the ECO II+ Tandem must be used, if necessary with a sufficiently dimensioned expansion container. The use of ECO II+ D / ECO II+ Tandem D is mandatory in Germany (see Appendix: General technical approval / DIBT).

26 Installation directly into the waste water pipe of the dental units

The ECO II+ is installed after the air/water separation (1.1) directly into the waste water pipe of the dental unit.

27 Installation after a central suction system (waste water < 1 l/min)

The ECO II+ is installed in the waste water pipe after the central suction system with water/air separation (waste water outlet).

28 Installation after a central suction system (waste water > 1 l/min)

The ECO II+ D (Germany) is installed in the waste water pipe after the central suction system with water/air separation (waste water outlet).

29 Installation after a central suction system (waste water < 2 l/min)

The ECO II+ Tandem is installed in the waste water pipe after the central suction system with water/air separation (waste water outlet).

30 Installation after a central suction system (waste water > 2 l/min)

The ECO II+ Tandem D (Germany) is installed in the waste water pipe after the central suction system with water/air separation (waste water outlet).

Installation in the suction line before a liquid ring pump

In liquid ring systems, the ECO II+ and the ECO II+ Tandem are installed in the suction line before the liquid ring pump.

3.3 Installation

31 Mounting module 1 (mounting element/wall bracket): Slide the two mounting elements together as in the picture.. The clips snap into place by pressing them together.

32 For the mounting the following minimal space is required (H x W x D): 360 x 165 x 202 mm. Position mounting element on an even floor space or fix onto wall with suitable screws.

33 Mounting module 2: Place module 2 on the device (36-1) and fasten it with the tensioning straps on both sides..

Hose connections:
34 **34-1** - Waste water inlet: waste water from the dental unit or from the central suction system;
34-2 - Waste water outlet: cleaned waste water directly to the sewage pipe

Installation ECO II+ versions:

35	36	37	38
-----------	-----------	-----------	-----------

 The waste water must be able to drain off freely into the sewage pipe! Secure fittings and tubes with clips!

3.4 Hose connections

39	40
-----------	-----------

A expansion container	D waste water connection
B ventilation valve	E suction system
C amalgam separator	F suction flow inlet

3.5 EXCOM hybrid suction systems with ECO II+ amalgam separation

- 41** EXCOM hybrid A1 - ECO II+ | EXCOM hybrid A2 - ECO II+
- 42** EXCOM hybrid A2 - ECO II+ D (Germany)
- 43** EXCOM hybrid A5 - ECO II+ Tandem D (Germany)

3.6 Filter

 A pre-filter with a mesh size of ≤ 4 mm must be installed upstream of the ECO II+. With METASYS suction devices, this is already permanently installed on the suction machine. The pre-filter must be cleaned daily.

4. Care and maintenance

4.1 Regular cleaning measures

i **After every treatment:**
 To remove residues from the lines and the suction system, the rinsing basin must be operated briefly after each treatment and each suction hose must be flushed with cold water.

GREEN&CLEAN M2:

2 x daily (noon/evening) and after surgical interventions, disinfection must be carried out with the prescribed disinfectant and cleaning agent GREEN&CLEAN M2. Ideally, GREEN&CLEAN M2 should be used prior to the treatment unit remaining at a standstill for extended periods (lunchtime, end of the day or holidays). For information on use and safety instructions, see the GREEN&CLEAN M2 instructions for use.

 It is not permitted to operate the ECO II+ amalgam separator without a pre-filter in the treatment units and/or the central suction system!
 The residues containing amalgam from the pre-filter can be collected in the METASYS ECOCENTER and properly disposed of with ECOTRANSFORM.
 The installed pre-filter must be cleaned daily.

4.2 Replacement intervals

The maximum service life (replacement interval) of the ECO II+ amalgam separator results from the chosen installation option and the number of operators.

Separator ECO II+	No. of operators	Max. service life (months)
1	1	12
	2	6
	3	9
2	4	6
	5	5
	6	4

 The resulting exchange interval must be observed! The ECO II+ must be replaced as soon as the maximum filling level has been reached, or 1x a year even if the maximum filling level has not yet been reached.

4.3 Filling level

44 A label on the ECO II+ collection container shows the filling level (100%, 95%) The label enables the user to monitor the filling level (see Appendix: General technical approval / DIBT).

 When using insoluble air abrasion powders the service life/replacement intervals need to be defined in the applicant's responsibility considering the expected consumption of the abrasion powder.
 The water level must not be exceeded during operation. If this mark is exceeded during operation, a service technician must be contacted.
 When the ECO II+ is not in operation, i.e. when the suction system is not running, the liquid level must be within the range of the „MIN. OPERATING WATER LEVEL“ mark. Check the operating water level daily.

4.4 Replacing the collection container

The ECO II+ collection container must be exchanged according to the replacement intervals as defined in the replacement interval table.

Risk of contamination:
 To avoid infection, wear personal protective equipment (hand, eye, nose and mouth protection) and disinfect and clean the device!
 Avoid contact with the contents of the amalgam separator!
 It is not possible to open the ECO II+ and reuse it! Every attempt to open the amalgam separator results in the malfunction of the device. Moreover, the required separation rate will no longer be achieved!

- 45** Remove the adapter with the tube at the waste water inlet and waste water outlet.
- 46** Tightly close the outlet of the full amalgam separator with the sealing plug (supplied with the amalgam separator) and secure with a clip.
- 47** Release the tensioning straps on both sides of the container (52-1) and remove the full amalgam separator upwards.
 Position the full amalgam separator horizontally onto an even surface, make sure that the inlet opening is pointing upwards!
- 48** Pour disinfectant for suction systems into the full collection container (20 ml GREEN&CLEAN M2).
- 49** Afterwards also tightly close the inlet opening with the sealing plug and secure with a clip.
- 50** Put the amalgam separator into the styrofoam shells, then place into the cardboard box.
- 51** Place the new amalgam separator on the fixture (56-1) and secure it with the tensioning straps on both sides.
- 52** Reconnect the adapter with the tube to the "IN" and "OUT" openings.

 Check inlet and outlet of the ECO II+ for tightness!

4.5 Disposal of the full collection container

Risk of contamination:
To avoid infection, wear personal protective equipment (hand, eye, nose and mouth protection) and disinfect and clean the device!
Avoid contact with the contents of the amalgam separator!

Certificates of recycling need to be stored according to national law.

Disposal with ECOTRANSFORM

53 Place the properly sealed amalgam separator into the two styrofoam shells.

Put the amalgam separator in the Styrofoam shells into the transport box (with postal certification) and close it.

54 Place the appropriate return label onto the transport box for recycling at METASYS logistics & collection GmbH (more information on the recycling options in your country on www.metasys.com).

5. Order numbers

101000022	Amalgam Separator ECO II+ international, intro kit, central
101000023	Amalgam Separator ECO II+ Germany, intro kit with expansion container, central
101000024	Amalgam Separator ECO II+ Tandem Germany, intro kit with expansion container, central
104000012	EXCOM hybrid A2 - ECO II+ 2 units, 230V, with amalgam separator ECO II+, 0.75 kW, 180 mbar
104000014	EXCOM hybrid A5 - ECO II+ Tandem 3-5 units, 230V, with ECO II+ Tandem, 1.1 kW, 180 mbar
104000015	EXCOM hybrid A5 - ECO II+ Tandem 3-5 units, 400V, with ECO II+ Tandem, 1.1 kW, 180 mbar
104000013	EXCOM hybrid A2 - ECO II+ Germany, 2 units, 230V, with ECO II+ D, with expansion container
104000016	EXCOM hybrid A5 - ECO II+ Tandem Germany, 3-5 units, 230V, with ECO II+ Tandem Germany, with expansion container
104000017	EXCOM hybrid A5 - ECO II+ Tandem Germany, 3-5 units, 400V, with ECO II+ Tandem Germany, with expansion container
104000011	EXCOM hybrid A1 - ECO II+ 1-2 units, 230V, 50/60 Hz, 0.75 / 0.9 kW, 180 mbar
120000630	Cover for ECO II+ / ECO Air incl. adapters
120000631	Expansion container for ECO II+ Germany / ECO II+ Tandem Germany

6. Warranty conditions

METASYS grants a guarantee of 12-36 months for specific products (duration of guarantee depends on the product according to the information in the current price list).

The guarantee covers all material faults that more than negligibly affect the function of the device. The guarantee does not cover damages caused by incorrect or improper handling as well as normal wear. Furthermore, the guarantee does not apply to the replacement of the amalgam collection container or to fragile parts such as glass, plastic, hoses, filters, condensate filters or membranes. Any incurred working and travel times are excluded from the warranty.

In order to determine the validity of the guarantee, the installation proof accompanying the device must be returned to METASYS immediately after proper assembly has taken place. In this case the warranty period begins with commissioning. In the event of installation without returning the installation proof to METASYS, any warranty claim is forfeited. Installation and return of installation proof must be completed within 24 months from the date of sale from METASYS.

Furthermore, any warranty claims of the customer expire if only one of the following circumstances arises, regardless of whether the circumstances arise for the customer of METASYS or a later owner or operator:

- Improper installation, operation, maintenance or transport of the device. If METASYS parts have to be returned, the original METASYS packaging must be used for shipping. Prior to packaging and shipping, the METASYS device must be cleaned and disinfected. Any openings where residual fluids could leak must be closed.
- Installation and return of installation proof have not been completed in the aforementioned 24 months.
- Failure to send the Installation Proof to METASYS.
- Installation and use of non-original METASYS parts.
- Installation of the device by personnel who are neither trained nor authorized by METASYS.
- Occurrence of damages through improper handling and operation or use of unapproved cleaning and/or disinfecting material, as well as non-compliance with the instructions for use.
- Execution of repairs by unauthorised repair shops or unauthorised personnel.
- Failure to comply with the prescribed maintenance intervals. Maintenance must be carried out 11-12

/ 23-24 / 35-36 months after installation of the respective METASYS item.

- Missing entries in the Equipment Logbook regarding the installation as well as prescribed servicing by technicians trained by METASYS.
- Failure to take reasonable immediate measures to avoid further damages in the event of a malfunction.
- Shipping of devices or components to METASYS without proper accompanying paperwork, in particular without error description or invoice for the purchase of the device.
- Failure to provide visual images (photos, video clips ...) of the METASYS item complained about, as well as of its installation situation and ambient conditions.

METASYS reserves the right to demand the documentation supplied with the device to check the maintenance intervals for the assertion of warranty claims. The processing of warranty claims takes place exclusively according to the following method:

In the event of malfunctions, the device must be opened by an authorized technician, the relevant component taken out and sent to METASYS unopened and cleaned. The customer sends the device or component in question to METASYS at his own expense. METASYS checks whether there is a warranty claim. METASYS will repair the device or component if it is cost-effective. The customer is charged for the costs incurred for the repair, but not the replacement parts covered by the guarantee. The consignment of the device or component to METASYS always represents a repair order for METASYS. For cost estimates for the repair of returned devices, a processing fee* is charged if the warranty period has expired or no warranty case exists. For products inspections only without any cost estimates a processing fee* can be charged. When sending the device or component to METASYS, an error description with all importation information about the device must always be included. The customer of METASYS may only provide payment in advance after consultation with METASYS. Only the affected component is to be sent in (smallest possible unit). If contaminated, intact parts are sent to METASYS without technical necessity, METASYS is entitled to destroy them without separate payment. The new part corresponding to the part to be destroyed is only delivered after a separate order and invoice. In any case, METASYS has the right to handle the guarantee through credit note or the return of new parts without conducting repairs. Guarantee services do not result in an extension of the warranty period nor do they initiate a new warranty period. The warranty period for installed replacement parts ends with that for the originally delivered device. The METASYS customer is obligated to make his customer aware of the conditions regarding the warranty processing. The statutory warranty rights of the customer remain unaffected.

* The current warranty conditions and fees can be found in the current METASYS price list.

Subject to printing and setting errors!

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